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Instruction for use of the device

Version 1.3—December 2025

This document is Mist Studio's property.

It may not be reproduced or communicated to third parties without the authorization of the top management.

Legende

Warning, consult associated documentation – **important safety rules**

Generalities and advice

Explanations

Population concerned

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Glossaire

Terme	Définition
Auroria game	Therapeutic game in virtual reality composed by a set
Auroria app	Mobile application allows to control remotely Auroria game. It allows also to display the content of the hel-
Level of Auroria game	A game level corresponds to a mission in the visual and auditory universe determined by the level. This is composed by different gate and a final enemy.
Intense immersive event	Content allowing higher immersion for the patient than the story mode. This is a final enemy or mini-game involving patient's cognition (reflex, mental arithmetic...).
Cast	Term used to define the action to display in live, the helmet contents. Here, the content is displayed to the tablet through <i>Auroria app</i> .
Virtual Reality (VR)	Virtual reality is an immersive technology that simulates 3D environments, allowing realistic interaction in a digital world.

Labelling of the product.....	42
Intended use.....	43
1—General Information.....	44-45
2— Safety instructions.....	46-49
3—Device description	50-51
4—Product installation	52-54
5—Stop AURORIA	55
6— Users.....	56-57
7—Patient account.....	58
8— Auroria game.....	59
9—Auroria app.....	60-69
10—Using accessories.....	70-71
11—Internet Connexion.....	72-73
12- Troubleshooting and ASV.....	74-75

Meaning of symbols and labeling

	AURORIA—Medical device class I
	<u>Manufacturer</u> : Mist Studio <u>Location</u> : 1 allé Saint Cloud, Villers-lès-Nancy, 54600, FRANCE <u>Legal representative of this company</u> : Matthieu Burtin, CEO
	Dispositif Medical certification
	CE certification
	Application version in I.X.X format Application Controller AURORIA : 1.0.X VR AURORIA Game : 1.0.X
	03770037649004AURORIASA
	Date of release (09-2025)
 eIFU indicator	Electronic instruction for use (e IFU)
	Indicates a website where a patient can obtain further information about the medical product. https://www.auroria.fr/

Intended used

AURORIA is a therapeutic virtual reality game (Auroria game) dedicated to the management of pain, soothing difficult care and, more generally, improving the care experience of patients in their medical pathway. Engaging in play helps patients divert their attention during medical procedures by focusing on game tasks, thereby reducing their perception of pain.

The product is specially designed to please pediatric patients, but AURORIA can be used safely for every patient aged more than 5-years-old.

The game is played by patients and controlled (real time content adjustment) by a caregiver (having completed training in the use of AURORIA) thanks to a mobile app (Auroria app).

AURORIA is used by painful patient during healthcare, at every stage of the patient's journey, in the operating room during surgery under local-regional anesthesia, in pre- and post-operative nursing care, during invasive examinations, or in the rehabilitation center.

WARNING !

Respect the restrictions and instructions provided in certain rooms/services of the establishment aimed at limiting or avoiding the use of radioelectric devices.

WARNING !

Respect contraindication

1.1 Software Configuration and Installation

The AURORIA software is installed on the proposed media. No configuration is required to launch the product.

1.2 Accessory Compatibility

The AURORIA software is compatible with the following accessories:

Virtual reality headset compatible with	Tablet compatible with AURORIA
- Minimum Chip: Qualcomm Snapdragon XR2 Gen 2 - Minimum Ram: 8GB - Minimum Storage: 128GB - Meta system - Minimum refresh rate 90 Hz - 1832x1920 resolution per eye - Maximum weight of 515g	- Minimum Chip or equivalent: Octa-core (1x3.36 GHz Cortex-X3 & 2x2.8 GHz Cortex-A715 & 2x2.8 GHz Cortex-A710 & 3x2.0 GHz Cortex-A510) - Minimum GPU or equivalent : Adreno 740 - Minimum Ram : 8GB - Minimum OS : Android 12

1.3 Performance du dispositif

- Pain reduction of up to 72% (Liu et al., 2021)
- 40% reduction in pain for the most painful medical procedures (Hoffman et al., 2019)
- Reduction in care time of up to 50% for pediatric procedural care (Jeff et al., 2024)
- In 65% of cases, the use of virtual reality games reduces the number of caregivers needed to carry out a complicated pediatric procedure (Chan et al., 2019)
- Caregiver satisfaction estimated at 90% (Le May et al., 2020)
- Parent satisfaction estimated at 90% (Savas et al., 2024)
- Patient satisfaction expressed by 94% of children in emergency wards who wish to use virtual reality games again during their future medical procedures (Chan et al., 2019)
- Calmer and less anxious patients with twice the fun of conventional care for children aged 7 to 17 (Atzori et al., 2022)

1.4 Lifespan of use

Accessories have a lifespan estimated to 3 years. Software and accessories are reusable.

i If a technical issue with the hardware or accessories occurs before then, please return it to our after-sales service. Any damage outside of the warranty period will be the customer's responsibility.

1.5 Environment of use

The product can be used in hospitals, clinical, emergency room or rehabilitation center.

WARNING !

Respect the restrictions and instructions provided in certain rooms/services of the establishment aimed at limiting or avoiding the use of radioelectric devices.

A stable internet connection is necessary to use AURORIA thanks to a WIFI network.

i Internet connection

It is necessary to connect the medical device to the internet at least once every 10 days.

1.6 Storage and preservation conditions

Tablet	VR Helmet
Water and dust resistant (IP68 rating)	Do not put in contact with water and liquids Do not expose the device to high concentrations of dust

i For more details, please refer to the user manual for this equipment.

1.7 Available language

The AURORIA product is available in French.

For more details regarding the operation of the associated equipment, please refer to the user manual for its equipment.

2.1 Security and maintenance

- Connect the devices (headset and tablet) to the internet at least every 10 days.
- Check for tablet and headset updates every 10 days.
- If a technical incident, a pairing problem between the headset and the tablet, or a connection problem occurs, press the "besoin d'aide" button at the bottom of the connection page.

2.2 IT-security measures

- The tablet must be protected with a password.
- AURORIA allows the creation of a user account for patient in order to have a follow-up concerning her/him statistics of game. No medical data are used or stored. The use of user account is not mandatory. User account will be protected by a unique identifier associated with a single patient.

2.3 Exclusion

AURORIA uses a virtual reality headset, which is a radioelectric accessory (a device that emits waves). The headset contains magnets and components emanating from magnetic/electromagnetic fields that may interfere with the operation of electronic devices and medical devices located nearby, including pacemakers, hearing aids and defibrillators. The product cannot be used during radiology examinations, MRIs, or in the presence of any other medical device that could interfere with the presence of the headset.

WARNING !

Observe the restrictions and instructions provided in certain rooms/services of the establishment aimed at limiting or avoiding the use of radioelectric devices.

2.4 Exclusion

 **WARNING** Patients with a history of epileptic seizures should not use AURORIA.

2.5 Contraindication

 WARNING !

Contraindications are as follows:

- Major psychiatric or neurological disorders
- Dizziness, or migraine headaches precipitated by visual auras including epilepsy
- Any wounds that may interfere with wearing helmet
- Visual, hearing or cognitive/motor impairments
- A significant reluctance to use the device
- Be less than 5 years old

2.6 Undesirable side effects

The residual risks claim as “side effects” are:

- Disorientation
- Dizziness, Motion sickness
- Fatigue, eye strain
- Headache
- Nausea
- Vision disorder/ other visual abnormalities
- Seizure

Risk control measures were implemented during product development. The benefit-risk ratio is acceptable for the product.

 Side effects

If a patient experiences a serious adverse effect, we invite you to inform the competent authority in your country and Mist Studio by email at support@auroria.fr

2.7 Precaution for use

Precautions for using the device are following.

Concerning helmet & controllers:

- Check the condition of the headset before using the product. If the headset appears damaged, do not use it.
- Do not bring the medical device into contact with water or liquids.
- It recommends that the helmet and tablet batterie should be completely charged before use.
- The product has to be disinfected before and after each use
- Do not charge or use your headset or controllers if you suspect or know that they are damaged.
- The helmet is a radio equipment. Keep away from the stomach of pregnant women and the abdomen of teenagers. Children and teenagers 7 should use radio equipment with moderation
- The helmet contains magnets and components emitting magnetic/ electromagnetic fields that may interfere with the operation of nearby electronic and medical devices, including pacemakers, 8 hearing aids and defibrillators.
- If the patient is pregnant, elderly, suffering from abnormal binocular vision, psychiatric or cardiac disorders or any other significant health problem, medical advice is required before using the medical device
- Do not charge or use your headset or controllers if you suspect or know that they are damaged.

Concerning the use of the product:

- The caregiver must position the headset so that the patient is comfortable throughout the treatment. The headset must not bother or hurt him. Special attention should be paid to children or people with muscle weakness regarding the weight of the headset.

- Always **use the medical device under the supervision** of a healthcare professional.
- The patient should not stand up with the headset on his head.
- Regularly check the temperature of your helmet and controllers to avoid burns and discomfort. Check mainly the part in contact with the patient's skin. The front of the helmet can potentially heat up without affecting patient comfort.
- During use, the caregiver must ensure that the helmet is tolerated by the patient.
- Always check the patient's pain. If the product appears to provide insufficient pain reduction, it is always possible to supplement the device with other analgesic techniques.

2.8 Hygiene

The device should be cleaned and disinfected by damp wiping the surfaces with a single-use disinfectant detergent wipe.

Warning

The device has electrical components: the wipe (or cloth) must be damp, not soaked. The wipe should be soaked, not the surface to be treated.

The device cannot be cleaned with a steam cleaner. UV cleaning is possible.

 Follow the establishment's procedure for cleaning non-sterile devices

3.1 Description of the product

AURORIA is a **medical device class I** aims to **reduce patient pain and to calm** the patient during painful or unpleasant care. It is an immersive therapeutic game that is used as a tool for health professionals to accompany patients during painful or unpleasant procedural care thanks to the immersive universe and adaptative gameplay.

This device is composed of two parts: a **virtual therapeutic game (Auroria game)** supported by a virtual reality headset


and a **mobile application (Auroria app)** available on tablet, aims to control the game.

This interactive game helps distract patients so that they feel less pain and therefore feel more appeasement. It is adaptive for patients thanks to an event system ("**AURORIA Focus**") controlled by the practitioner allowing them to overcome a potential acute pain.

3.2 Accessoires

Check the accessories and documentation following this list.*

VR Helmet	Controller, battery and protection	Associated power cable	Silicone Protection	Tablet	Tablet shell
					
External	Strap of the helmet				
					

Documentation :

Instruction for use	Memo	Tutorial Video	Instruction for use—helmet	Instruction for use—tablet
				

*The images presented in this notice are for illustrative purposes and may differ from the version provided.

i Note: hygiene

Before and after each use of the medical device, it is necessary to clean the devices in order to allow their disinfection.

i Follow the installation steps in the order given

4.1 Preparation of the tablet

- 1) **Unlock** the tablet by pressing button (A)
- 2) Write **password**
- 3) Check **internet connection** (refer to p.30)
- 4) Open **AURORIA application** (B)

5) Chose the **wanted mode**. We recommend using guest mode if the patient is not recurring. .

An **adress IP** is asked. To find it, prepare the virtual reality helmet. Instructions are given in the part « helmet » of the application

i Internet network: the network must be stable, and must be the same between the tablet and the headset.

4.2 Preparation of the casque


Turn on the headset

- 1) To switch on the headset, long press the button located on the left part, on the front of the helmet.
- 2) Manually define the play area ("**Roomscale**") by drawing an area that includes the patient and caregiver.

3) (**OPTIONAL**) Enable mirroring (display the headset's content on the tablet)

1. To do this, tap Settings (bottom left of the screen)
2. Select mirroring
3. Choose "Web"
4. Click "Next"
5. A code will appear; enter it on the left side of the tablet

4) Turn on **Auroria application** by pressing Auroria application.

At this moment, an IP adress is displayed. It must be entered on the mobile application AURORIA on tablet.

i Follow the installation steps in the order given

4.3 Preparation of the patient

- 1) Put the patient in the **position of the care, head straight**
- 2) **Recenter the position** (align patient vision and game content) with the button home (G) of the right controller by pressing it during 5 seconds.
- 3) **Loosen** the helmet strap as far as possible by unscrewing the knob (H) on the back of the helmet.

- 4) Put the helmet on the patient's head, starting with the eyes. Depending on the frame, it may be possible to keep the glasses on during play.

- 5) Tighten the knob (H) behind the head until the helmet fits snugly.
- 6) Place the associated controller in the patient's hand (if right-handed, give the right lever; if left-handed, give the left controller).

⚠ The patient should not move while wearing the helmet.

⚠ Always ensure that the patient is comfortably installed

5.1 Stop Auroria game

5.1.1 Stop current content

Retour au QG

To stop the game, you can select "**Retour au QG**." The patient will then be sent to a **transition zone without** any **distracting** events. This corresponds to the player's starting zone (refer to the "Game Level" document for more information about QG).

5.1.2 Stop use AURORIA

Terminer le soin

To stop using the AURORIA product, including the Auroria **game**, press the "Terminer le soin" button on the "Menu" page of the Auroria app.

This gently **shuts down the headset** so the patient can return to the real world. Remove the headset from the patient once prompted to do so on the mobile app.

i End of the game

A few second is necessary between the trigger of the end and possibility to remove the helmet. This avoids a potential disorientation for the patient due to a sudden interruption of the virtual reality

5.2 Remove the helmet

Afterwards, you can remove the headset from the patient. To do this:

- 1) Retrieve the patient's game controller.
- 2) Loosen the wheel on the virtual reality headset and remove the headset from the patient.
- 3) Turning off the headset is done by pressing the button on the left side of the front of the headset.

i After each using

- Clean the equipment with the usual cleaning tools.
- Charge the headset and tablet.

6.1 Intended users

The video game **Auroria game** is used by **patients** during a care, and **controlled** (choice of level, immersive event triggering) by **caregiver** thanks to the mobile application associated **Auroria app**.

Accompanying person of patient can use the application's cooperation mode to play a game and send a gift to the patient in the game.

Caregiver	Patient	Accompanying person of patient
Nurse, doctor, surgeon who have validated training for the use of AURORIA	Patients without contraindication who have undergone a painful or unpleasant medical procedure and are over 5 years of age.	Parents or accompanying must be able to see and read to use the product. There is no special knowledge or skill required to use AURORIA No training is required to use the video game part of the application. They can only use this part of the application.

 **AURORIA** is written in French. Caregivers have to be able to read and understand indications displayed on the tablet screen, and patients have to be understood French oral instruction in **Auroria game**. No prerequisites concerning skills or knowledge are necessary for **AURORIA**.

 **WARNING !** A training is required to use the product Auroria app.

The training, supervised by a member of Mist Studio team, is performed by each caregiver that will use AURORIA. This training ensures a good use of the product.

Training in the use of AURORIA is **not associated with continuing vocational training** and is not included in the skills development plan. The purpose of the training is to teach how to use the product.

6.2 Characterization of the patient population

Every patient aged more than 5 years old can be used **AURORIA** to reduce pain and anxiety, before, during and after a medical procedure.

Attention to the patient's mental and physical abilities

Caregivers in charge of the intervention can ensure the patient's mental and physical abilities are appropriate for the use of the product. Patients with psychological, psychiatric, cognitive, visual, hearing disorders and photosensitive epilepsy, should refer to the medical advice of the physician in charge of the care or/and the patient generalist practitioner before using the **AURORIA**.

No prerequisites regarding the use of virtual reality, knowledge or skills related to games are necessary for the use of the product, although an appetite for gaming is necessary. If the patient refuses the use of **AURORIA**, carry out the classic procedure.

Attention to the patient's injuries


If the patient has lesions in the area above the neck, refer to the medical advice of the doctor in charge of the procedure and/or the treating physician for the use of **AURORIA**.

The product can be used when the patient is sitting, lying on his stomach or in a semi-recumbent position. Adjust the patient's position so that it is comfortable.

7.1 Create an account

Each patient can have a user account.

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7.2 Connect to a patient account

If patient already has an account, write the patient ID generated during the creation of the account (Page « Nouvel identifiant »).

7.3 Invited mode

It is advised to use “invited mode” if the product will only be used once by the patient. In this case, no follow-up and progression of the patient play will be registered.

It is advisable to note the identifier generated for the patient.

No information regarding the patient's identity is stored by the product.

8.1 Presentation

?	Auroria game is the therapeutic game used by the patient to reduce their pain and anxiety during medical care.
	Used by patients aged over 5 years who have no contraindications
	1. Shoot enemies to fight them 2. Shoot targets to earn bonuses and points

8.2 Instruction of game

The game instructions are given directly in the content sent to the headset. The patient will therefore be able to understand what is expected of him thanks to the game audio.

Help from the caregiver

At the same time, the caregiver has access to instructions or help to give to the patient if he feels the need. (see p.24)

8.3 Game Control

Shooting an Enemy


To shoot an enemy, simply press the round button (A) on the back of the controller.

9.1 Presentation

	Auroria app is the mobile application available on tablet which controls Auroria game .. It serves as a sort of "remote control" for the game.
	Used by caregivers and parents or companions of the patient.
	1. Show headset content 2. Control the elements sent to the headset 3. Film the caregiver/accompanying person to reassure the patient

9.2 Description of the interface

« helmet » part

allowing the contents of the headset to be displayed

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allowing content to be controlled

The tablet screen, once the **Auroria app** is launched, is divided into **two parts**.

The right part corresponds to the “remote control” part of the application. This part allows the patient to connect to their account (1), and is intended to control the content sent to the headset (2).

The **left part** corresponds to the contents of the headset, available either by:

Mode direct

allowing the headset content to be **displayed live**, thus allowing the patient's interactions to be seen. This mode requires a stable internet connection.

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Recommended mode for displaying the player's progress through the content they are exploring as **images**. Information about the world and the tasks to be completed will be displayed.


“Direct Mode” requires upstream mirroring. For this, refer to step 3) p.13 “Headset Installation”.

If mirroring is enabled, it is possible to switch between modes while using the product.

9.3 Buttons

Fonction-	Fonctionnelles	Remarque
Launch a content		
End of care	 	Completely stop care: the game is gradually cut to avoid patient disorientation
Trigger an emergency event	 « Warning button »	Located at the bottom of the screen: available at every moment, to choose between immersive event or camera activation.
Launch an intense immersive event		Accessible thanks to the « Warning button ». This button triggers randomly an intense immersive event.
		Thanks to home page, it is possible to choose intense immersion event which be send to the patient.
Change the game control mode		Allows you to choose whether the player aims at targets with their head, the left controller, or the right controller (which makes a weapon appear in their hand).
Cooperation Mode		Access to the CO-OP page
Back to the QG		Allows the patient to return to a transition area, without any particular distractions.
Help		On the login page, this button provides access to the online help center (user manual, video, customer service).

9.4. Different modes of game

9.3.1 The game control mode


Aiming the weapon allows the player to be more immersed, but they must be able to move their hand and arm during the treatment. The weapon represented in the game corresponds to the controller that the patient holds in their hand.

Aiming with the head allows the patient to keep their arm static. Aiming is achieved solely by head movement. To shoot, simply press the trigger (see 10.1 Using the head-set) without moving your arm or hand.

9.4. Differents mode of game

On the home page, it is possible to choose between 3 options:


	Mode Aventure	Mode Défi	Liste des Focus
What ?	Mode par défaut Allows you to select	story mode level revisited	Distracting events categorized
Who ?	All patients	Experienced patients	All patients
When ?	From the beginning of the treatment	From the beginning of the treatment	In anticipation of a painful gesture or to calm an anxious state

9.5 Description of a content page

Summary map of the content.


If we click on it, it takes us to the presentation of this content

The different types of content are explained on page 23

9.6 Differentes pages of the application

9.6.1 Home page

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
Mode Aventure

The **Mode Aventure** correspond at all levels of AURORIA covering the different universes.

The list of game content is present in the attached document « AURORIA contenu »

Return button

Card corresponding to a content
(Displays content description if we click on)

All pages are structured in this way


9.6.3 « AURORIA Focus» Page

Liste des Focus

« **Liste des Focus** » gather the game's most distractive events together. This allows caregivers to intervene when patients are too agitated or in pain to proposing a game level more intense in terms of immersion. This is possible thanks to contents call on cognition and/or which ask a lot of interaction with the virtual word, for example. Moreover, a change of rhythm during the game « **Mode Aventure** » avoids patients dropping out, distracted by the new virtual environment around them.

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Mode Défi

The « **Mode Défis** » page retakes the same level of the "Mode Aventure" page (with higher level of difficulty. This allows to the most experienced gamers, or patients who want more difficulties, to access to a richer gameplay, keeping AURORIA universe.

9.6.5 « Co-op» page

The CO-OP part allows patient companions to play games during patient care. This way, the companion can send live bonuses to help the patient achieve their gaming goals.

Always ensure that the patient is comfortably installed

9.6.6 Description of the page « Co-op»

Interface of the game « Co-op »

Amount of points collected


Once the yellow gauge is filled, a bonus level is unlocked

Note: When the game ends, the number of points is kept.

Enemy: Touch him to gain points

Golden enemy: gains more points but only stays 3s

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Retour

(1) Shoot enemies

Press them directly to collect points

(2) Redeem your points

Use your points by pressing the button located at the bottom right of the screen.

(3) Unlock a Bonus Level

Reach the goal in the top left corner of the screen to unlock a bonus level in the shop menu.


(4) Stop the game

Quit the game at any time by tapping the bottom left corner.

10.1 Helmet using

Action	Description	Related Illustrations
Turn on or turn off the helmet	To turn on the headset, long press the button on the left side, at the front of the headset (1)	
Display the menu	If you have any problems, you can display the headset menu by pressing the Meta button (2) on the right controller.	
Validate an action	To validate an action on Meta, press the trigger located behind the controller (3)	
Calibrate the screen	If the image displayed on the screen is not centered in relation to the user, hold down the Meta button (2) on the right controller for 5 seconds.	

10.2 Tablet using

Action	Description	Related Illustrations
Turn on	Long press the top left button (button A) of the tablet	
Turn off	Slide your finger on the right of the screen, from the top to the bottom to display the menu Press the power off button at the top of the screen (C)	

i For more information, refer to the accessories user manual.

Internet connection for using AURORIA is necessary to pair the both device (tablet and helmet).

i The tablet and the helmet have to be connected on the same network.

11.1 Access to WIFI information and connexion


To check connection regarding the tablet:

Check that the wifi logo is white and not gray (A)

Check the name of the network it is connected to (displayed below the word "Wifi" (A.1))

To check connection regarding the helmet


- Go to headset settings
- (J) displays Wifi status (connected or not, or which network it's connected to)
- Press this button if you're not connected to the right network

11.2 Connection by SIM card

Mist Studio offers a tablet with the possibility of receiving a SIM card **allowing access to the internet from the tablet** without being connected to wifi.

- Check if the WIFI is inactivated: the (A) logo must be grey
- Check if internet connection is activated: the **logo (B) must be white** (means active).
- Activate mobile data allowing share connection to the helmet: the (C) logo must be white.

11.3 Wifi Connection

Refer to the person in charge of the institution.

Check before using AURORIA that the tablet and the helmet is connected on the same network. (A) logo must be white.

12.1 Concerning AURORIA

If the patient has no affinity with the world of video games, the product cannot be used. In this case, the classic care procedure must be followed.

Adverse effects

If a patient is subject to a serious adverse effect, we invite you to inform the competent authority of your country and the company Mist Studio, by email at support@auroria.fr

12.2 Concerning the helmet – Auroria game

The headset does not start	Check that the headset is sufficiently charged. If this does not work, contact customer service
The Auroria game does not start	Restart the device Check internet connections If it doesn't work, contact customer service
I can't find the app Auroria game	Check that the application is still installed by going to the tablet's search bar and typing AURORIA »
I can't connect to mobile data	Refer to after-sales service

12.3 Concerning the tablet – Auroria app

The tablet does not start	Check that the tablet is sufficiently charged. If this does not work, contact customer service
The application Auroria app doesn't start	Restart the device Check internet connections If it doesn't work, contact customer service
I can't find the app Auroria app	Check that the application is still installed by going to the tablet's search bar and typing « AURORIA »
I can't log in to the patient account	Create a new account or choose guest mode
I can't connect to mobile data	Refer to after sale service

After sales service information

If the game **Auroria game** or the application **Auroria app** show bugs, moments of latency, or other anomalies, contact support@auroria.fr

After-Sales Service Shipping

Before sending the equipment to Customer Service, please gather all the items in the kit and place them in the storage bag provided.